

Provitro AG

Charitéplatz 1

tel +49.30.450 578 358

sales@provitro.de www.provitro.de

Charité Campus Mitte

10117 Berlin fax +49.30.450 578 919

hMSC proliferation medium, basal

Cat.-Nr.: 200 0901

contains of:

Basal media		Supplements
200 0901	500 ml hMSC proliferation medium, basal	-

Maintenance of hMSC proliferation medium:

Place the bottle of **basal medium** in the dark at **4°C to 8°C** immediately after delivery.

Characteristics:

The Provitro hMSC proliferation medium is a sterile liquid culture medium for culturing human mesenchymal stem cells (hMSC). The medium is delivered as a basal medium and is suitable for culturing human mesenchymal stem cells **after adding optional available essential supplement mix components.** The final formulation is optimized for initial seeding of 6,000 cells / cm² up to confluence (approx. 90 %). Feeder-layer, matrix substrates or other substances are not necessary.

Stability and storage:

The supplemented hMSC proliferation medium can be stored in the dark at 4°C to 8°C for up to 1 month. Do not heat the medium over 37°C or use uncontrollable sources of heat (e.g. microwave appliances). If only a part of the medium is to be used, remove this amount from the bottle and heat it.

Special note:

Do not freeze the medium. This can lead to high salt concentrations by freezing out pure water which will cause irreversible damage.

Quality control:

Provitro's hMSC proliferation medium is thoroughly tested after each production. All components are tested in a stringent biological assay. Each batch is checked for hMSC proliferating characteristics. The cells cultured in hSMC proliferation medium are checked regarding their morphology, the adherence rate, the colony forming efficiency and the population doubling time.

Product specification:

The pH is set at 7.6 and osmolality at $285 \pm 10 \text{ mOsm} / \text{kg}$.

In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.