

Provitro AG Charité Campus Mitte Charitéplatz 1 10117 Berlin tel +49.30.450 578 358

fax +49.30.450 578 919

sales@provitro.de www.provitro.de

# Supplement kit for endothelial cell proliferation medium, HuS-kit

Cat.-Nr.: 216 0001

contains of:

Basal media	Supplements	
-	222 1000 L-Glutamine	
	232 3500 HuS (human serum AB)	
	226 0500 Heparin	
	244 0500 human rec. EGF (epidermal growth	າ factor)
	245 0250 human rec. bFGF (basic fibroblast)	growth factor)
	241 0025 human rec. VEGF	
	242 0500 human rec. Long R3 IGF-1	
	223 0005 Ascorbic acid	
	224 0010 Hydrocortisone	
	236 0350 Antibiotics (optional)	

#### Maintenance of supplement kit for endothelial cell growth medium:

The supplements are delivered on dry ice. Store the supplements at -20°C.

#### **Characteristics:**

The Provitro supplement kit for endothelial cell proliferation medium is suitable for culturing Provitro human endothelial cells after adding the supplement kit components to 500 ml basal medium. The formulation is optimized for initial seeding of 4,000 cells / cm² up to confluence (approx. 90 %). Feeder-layer, matrix substrates or other substances are not necessary. Due to the possibility of reduced proliferative activity we recommend to use the antibiotic supplement for freshly isolated cells only.

#### Stability and storage:

Store the **supplements** at **-20°C**. Do not use the supplements after expiring date.

### Special note:

Avoid any thaw/freeze cycle.

#### **Quality control:**

Provitro's supplement kit for endothelial cell priliferation medium is thoroughly tested after each production. All components are tested in a stringent biological assay. Each batch is checked for human endothelial cells proliferating characteristics. The cells cultured in endothelial cell growth medium are checked regarding their morphology, the adherence rate, the colony forming efficiency and the population doubling time.

## In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.