

Provitro AG

Charitéplatz 1

tel +49.30.450 578 358

sales@provitro.de www.provitro.de

Charité Campus Mitte

10117 Berlin fax +49.30.450 578 919

# Supplement kit for defined fibroblast maintenance medium, serum-free

Cat.-Nr.: 217 0401

contains of:

Basal media	Supplemen	Supplements	
-	222 1000	L-Glutamine	
	245 0050	human rec. bFGF (basic fibroblast growth factor)	
	246 0250	human rec. Insulin	
	236 0350	Antibiotics (optional)	

## Maintenance of supplement kit for defined fibroblast growth medium:

The supplements are delivered on dry ice. Store the supplements at -20°C.

#### **Characteristics:**

The Provitro supplement kit for defined fibroblast maintenance medium is suitable for maintenance culturing of Provitro human fibroblasts (HFIB) **after adding the supplement kit components to 500 ml basal medium**. The formulation is optimized for initial seeding of 6,000 cells / cm². Feeder-layer, matrix substrates or other substances are not necessary. **We recommend using the antibiotic supplement for freshly isolated cells only.** 

## Stability and storage:

Store the supplements at -20°C. Do not use the supplements after expiring date.

#### Special note:

Avoid any thaw/freeze cycle.

### **Quality control:**

Provitro's supplement kit for defined fibroblast maintenance medium is thoroughly tested after each production. All components are tested in a stringent biological assay. The cells cultured in defined fibroblast maintenance medium are checked regarding their morphology, the adherence rate, the colony forming efficiency and the population doubling time.

## In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.